

Adjusting a 'Select' ROHO® Cushion

1 Ensure the ISOFLO® Memory Control™ Unit valve is in the unlocked position by sliding the valve to the right.



2 Inflate the cushion until all the cells are full of air. When this is done the cushion will have a SLIGHT upward arch.



3 Place person on cushion. Check footplates are properly adjusted.



4 Slide your hand palm side down to locate the bony prominence—usually the ischial tuberosity. If you are not sure lift the persons leg slightly and you will feel the boney area move.



- Check daily for improper inflation, but adjust air only when required.
- Due to the design of the product very little, if any, air will leak unless there is a puncture.

Quadtro Select or Contour Select


5 While keeping your hand under the bony prominence, open the ROHO air valve and let air escape. Continue to let air out until you can feel about 15mm between the bone and the base of the cushion. The person will feel themselves sink into the cushion.



6 Close the air inlet valve. Rotate clockwise to close.



7 Position the pelvis. If correcting asymmetry you may need to hold the pelvis in position. Once the pelvis is in the desired position push the ISOFLO® Memory Control™ Unit to locked position.



'Contouring' (as per this photo) after a person gets off the cushion is normal.



If cells look full with nobody on the cushion it is likely it is over-inflated.

- Remember 'in' it - not 'on' it
- Placing things like sheepskin over the cushion will reduce the performance of the cushion.
- Do not use the ISOFLO switch or the cells as a handle. The yellow cord is there for this purpose.