

Universal Concave Cover For ProCair Range



- A valuable accessory for the ProCair, ProCair Plus, ProCair Trio, and ProCair Quad alternating pressure mattresses
- Helps prevent patients falling from the bed or injuring themselves from rolling against side rails
- Modular foam sections are easy to remove for cleaning, while the hinged sections conform to profiling beds
- Components are heat sealed to protect against moisture ingress for maximum infection control
- Universal fit with elastic corners and straps to secure the cover to an existing alternating pressure mattress

THERAPEUTIC DEVICES SHOULD ONLY BE USED IN ACCORDANCE WITH MANUFACTURERS INSTRUCTIONS AND UNDER THE CONSENT, SUPERVISION AND MANAGEMENT OF A SUITABLY QUALIFIED HEALTH PROFESSIONAL.

	BSCS1001	Universal Concave Cover, Single 2000 x 880mm
Code	BSCS1002	Universal Concave Cover, King Single 2000 x 1050mm
Warranty		2 years

Scan QR code for more information
novis.com.au/product/procair-concave-cover/

