

Over Bed Table

The BetterLiving Over Bed Table (OBT) is designed to provide a portable and stable surface that can be positioned over a bed, chair or similar. It features 4 locking swivel castors for easy maneuvering, spring loaded height adjustment and a thermoformed top for improved infection control as well as a raised edge to catch spills and help prevent items slipping off the table top.



NAME	BetterLiving Over Bed Table	
TOP COLOUR VARIATIONS	White	BL9990
	Off White Wood Effect	BL9991
	Brown Wood Effect	BL9992
OVERALL HEIGHT	715-1100 mm	
OVERALL WIDTH	810 mm	
OVERALL DEPTH	450 mm	
WIDTH BETWEEN BASE	672 mm	
TABLE SIZE	810x450 mm	
TOTAL WEIGHT	≈12 kg	
SAFE WORKING LOAD	20 kg	
WARRANTY	1 Year	


*measurements subject to change

Recessed Top




Collects and prevents liquids from spilling over the edge

Available Tabletop Finishes



White



Chocolate Brown Wood Effect



Off White Wood Effect