

7 Reasons to Consider The 'Casero Duo' Home Care Bed



1. ACCESSORIES/ AESTHETICS

Headboards, footboards, padded side panels, rails and side grip handles are all quoted with bed price as standard on the Casero Range.



2. OPTIONS

2 headboard designs and 4 fabric colour choices.



3. SWIVEL CASTORS

8 super strong swivel castors with large foot lock pedal.



4. ELECTRICS

Linak electrics with corded hand control. Ability to program a bed scanner with Sip n Puff, Joystick or Jelly Bean button.



5. BED LENGTH EXTENSION

Casero bed range has the ability to lengthen to provide an additional 100mm and 200mm. No additional pricing.



6. BED HEIGHT RANGE

Optimal home care bed height range 184 - 760 mm.



7. SINGLE WEIGHT CAPACITY

Single and King Single beds - 250 kg



The 'Casero Duo' Home Care Bed

Specifications

Dimension

Single

Patient Surface Width	: 900mm
Patient Surface Length	: 2000mm 2200mm
Length Extension	: 200mm
Overall Width (includes side surround)	: 1020mm
Overall Length (includes Head and Footboards)	: 2170mm 2370mm

King Single

Patient Surface Width	: 1050mm
Patient Surface Length	: 2000mm 2200mm
Length Extension	: 200mm
Overall Width (includes side surround)	: 1190mm
Overall Length (includes Head and Footboards)	: 2170mm 2370mm

Technical Specifications

Backrest Adjustability	: 0-65 Degrees
Knee Bend Range	: 0-32 Degrees
Trendelenburg Range	: +-12 Degrees
Height of the Base	: Adjustable from 184 to 760mm
Backrest Length	: 770mm
Fixed Section	: 280mm
Upper Leg Section	: 340mm
Lower Leg Section	: 660mm

Safety

Maximum Safe Working Load	: 250 kg
---------------------------	----------

Warranty

Frame	: 10 years
Electrics	: 5 years
Hand Control	: 1 year
Surrounds (includes Head and Footboards)	: 1 year

Compliance

Australian Standards	: IEC 60601-2-52
----------------------	------------------



Accessories



BZ0548 CASERO CLAMP ON RAIL (400MM)

Bed grab bar to help people get into and out of the beds. Installed securely on the bed frame and provide the user with a secure handhold.



BZ0443 CASERO TELESCOPIC SIDE RAIL (PAIR)

The European designed, vertically adjustable side rail options offer versatility, support, safety and protection. The bed can have no rails or any number 1-4 split rails.



PROCAIR RANGE ALTERNATING PRESSURE MATTRESSES

Novis offers a range of alternating pressure mattresses, hybrid foam air and foam mattresses.



01062023F