

# Invacare CS8 floorline bed

The Invacare CS8 bed is a fully-featured electric floorline bed with one button chair positioning, attendant panel, under bed lighting, backup battery and an ultra low 150mm low height.

## Features :

- Total safe working load (SWL) of the bed 250Kg.
- Linked braking on head and foot ends
- Electrically adjusted height range. The bed works within a height range of 150mm low height to a maximum high height of 680mm.
- Trendelenberg head and foot down tilt 14°.
- Four sectioned electronically operated mattress platform with independent adjustment of backrest angle (maximum 65°) and knee section angle (maximum 35°).
- Auto-regression backrest for reduced shear and friction.
- Integral foot end mattress base extension. Extends the bed by 195mm.
- (Optional) bolt on 3/4 length folding side rails.
- Removable timber head and foot panels.
- One button chair positioning.
- Under bed lighting.



# Invacare CS8 floorline bed

The Invacare CS8 bed offers a range of positional movements and additional features including auto-regression backrest and linked braking that assist caregivers and provide for patient comfort. The bed includes as standard: back-up battery, under bed lighting, and attendant control panel.

## Technical Specifications

<b>SIZE OPTIONS:</b>	Single	King Single
	W420-08-BB-V2	W420-08-KS-BB-V2
<b>WARRANTY:</b>	10 year frame   3 year electrical	
<b>BASIC DIMENSIONS:</b>		
OVERALL LENGTH	211.0cm	211.0cm
OVERALL LENGTH (BASE EXTENSION EXTENDED)	230.5cm	230.5cm
OVERALL WIDTH	92.0cm	109.0cm
OVERALL WIDTH (RAILS FITTED)	100.0cm	120.0cm
MATTRESS SURFACE LENGTH (BETWEEN PANELS)	199.5cm	199.5cm
BED PLATFORM HEIGHT RANGE	15.0cm to 68.0cm	15.0cm to 68.0cm
BACK REST ANGLE ADJUSTMENT	0° to 65°	0° to 65°
KNEE BEND ANGLE ADJUSTMENT	0° to 35°	0° to 35°
TRENDELENBERG TILT ADJUSTMENT (HEAD AND FOOT DOWN)	0° to 14°	0° to 14°
<b>MATTRESS SIZE:</b>		
RECOMMENDED MATTRESS SIZE	WIDTH	86.0cm
	LENGTH	198.0cm
MAXIMUM MATTRESS DEPTH (WHEN SIDE RAILS IN USE)	17.0cm	17.0cm
MINIMUM MATTRESS DEPTH (WHEN SIDE RAILS IN USE)	10.0cm	10.0cm
<b>SAFE WORKING LOAD:</b>		
S.W.L SAFETY FACTOR = 2 X S.W.L	250kg	
<b>ACCESSORIES (1 YEAR WARRANTY)</b>		
3/4 LENGTH FOLD DOWN SIDE RAILS, PAIR	CS8-FSR	
SELF HELP POLE	SP002	
IV POLE	IV110	
ASSIST BARS, PAIR	CS8-ASB	

