



Unique design of this seat evenly distributes body weight and encourages proper alignment of the pelvis and thighs for extended sitting comfort.

Use it with any ObusForme Backrest Support to create a complete and comfortable ergonomic sitting experience.

Made of soft, comfortable polyurethane foam with fabric cover, the seat relieves fatigue from the buttocks to the sacrum to the lumbar and thoracic areas.

- ◊ ideal for periods of prolonged sitting and to help prevent pressure injury

- ◊ polyurethane foam absorbs compression and dissipates vibration created in moving vehicles
- ◊ easily removable front section for use with narrow chairs and benches
- ◊ lightweight and portable for easy transport; includes elastic positioning strap for secure attachment

CODE	OBUSS
COLOUR	Black
WIDTH	450 mm
HEIGHT	400 mm
DEPTH	60 mm
MATERIALS	Polyurethane foam
WEIGHT	800 g
CAPACITY	100 kg
WARRANTY	1 Year



ENDORSED BY



Canadian
Physiotherapy
Association

