


Adaptor® Pad

All purpose,
cut to fit air cushion



Ideal for anywhere shock
protection is required

The self-contained AIR SPRING® ADAPTOR Pad can be used as an interface on virtually any flat or contoured surface to protect soft tissues. Especially useful when cut to fit various surface shapes such as sides of wheelchairs, trays, head and foot rests, inside of helmets, or anywhere uncontrolled movement may need protection.

- › cut to fit any flat or contoured surface
- › use as an interface to protect soft tissue, absorb shock and protect against shear and friction
- › ideal for head and foot rests
- › no cover, repair kit or pump
- › suitable for users of any weight

Code	R-ADAPTOR
Capacity	Unlimited
Width	340 mm
Depth	250 mm
Height	25 mm
Weight	1 kg
Conforms to	ARTG 111106
Warranty	2 year