

Universal Cordless Pads, Mats & Transmitter



FEATURES

Guaranteed **18 month warranty**

Universal Bed Pads are foldable (single fold)

Simple operation and maintenance

When fitted with a Universal Transmitter (2720), Pads & Mats work with cura¹ Cordless Receivers eg. 2700, 2544, 2549, 3500, 3501 & 3508

When fitted with a USB to RJ10 Breakaway Cable (2693), Pads & Mats work with a Premium Monitor or receiver eg. 2701, 2544

Easy to clean non-allergenic soft PVC

Floor Mat has low profile edges to reduce tripping hazards & allow easy wheelchair access

Foldable Floor Mats in beige colour to help some Alzheimer patients who have trouble with dark mats being perceived as a "dark hole in the floor"

TECHNICAL SPECIFICATIONS

Recommended weight range for reliable operation: 30 - 200 kg.

Pad / Mat Dimensions:

Chair Pad - 380 x 250 mm

Bed Pad - 760 x 250 mm

Floor Mat - 1220 x 610 mm

Floor Mat Single Fold - 900 x 600 mm

Floor Mat Double Fold - 1220 x 610 mm

Universal Transmitter Insert

Range: 50 m*

Operating frequency: 433 MHz

Incorporates status / in-use LED indication

ACCESSORIES

2700 Six Channel Cordless Monitor, **2701** Falls Monitor with SafeTime Support (corded connection only),

2720 Universal Cordless Transmitter, **3500** ActiveCare Console, **3501** ActiveCare Pager, **3508** LED Pager, **2693** Mini USB to RJ10 Breakaway Cable, **2544 & 2549** Multiport Expanders

cura¹.com

cura¹ products are designed to augment well designed care programs and practices, cura¹ does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision, and are at risk.

* Dependent on environmental conditions