



MASTERPIECE COLLECTION

ALVIO LIFT CHAIRS

Improving your
life everyday™

KCare
Healthcare Solutions

**The Alivio Masterpiece
Collection Lift
Chair range are a
combination of
style and support
with comfort**



PICASSO



LEONARDO



MICHAELANGELO



DONATELLO



**DONATELLO
PETITE**

PICASSO LIFT CHAIR

Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy to use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower the head position



Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight
Picasso Lift Chair KA551	1003mm	813mm	1080mm	508mm	521mm	495mm	1830mm	74kg	158kg

FOOTREST EXTENSION



LEONARDO LIFT CHAIR



Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy to use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower the head position

Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight
Leonardo Lift Chair KA552	1067mm	902mm	1092mm	521mm	571mm	483mm	1830mm	76kg	158kg

LAY-FLAT RECLINER



MICHELANGELO LIFT CHAIR



Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy to use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower the head position

Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight
Michelangelo Lift Chair KA553	1080mm	890mm	1067mm	546mm	559mm	470mm	1830mm	70kg	158kg

LUMBAR SUPPORT TO ACCOMMODATE SITTING POSTURE



DONATELLO LIFT CHAIR

The **Donatello** and **Donatello Petite** have a range of configurable back rest options to suit you.



**Lateral
Support
Backrest**
KA554-LS



**Standard
Support
Backrest**
KA554-SB



**High
Support
Backrest**
KA554-TB

DONATELLO PETITE LIFT CHAIR



**Lateral
Support
Backrest**
KA555-LS




**Standard
Support
Backrest**
KA555-SB

Features

- Two stage lift function
- Lay-flat chair
- Footrest extension
- Easy to use hand controller
- Ultra durable stylish fabric
- Lumbar support to accommodate sitting posture
- Power headrest to lift or lower the head position

Description	Overall Depth	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Reclined Length	Product Weight	Max User Weight
Donatello Lift Chair KA554	1016mm	890mm	1092mm	520mm	520mm	508mm	1830mm	74kg	158kg
Donatello Petite Lift Chair KA555	978mm	851mm	1041mm	457mm	457mm	457mm	1830mm	70kg	158kg

ULTRA DURABLE FABRICS



Four independent motors for a wide range of position combinations and support



Power Headrest

The individually adjustable head support allows you to lift or lower the position of the headrest even when the backrest is fully reclined providing a higher level of adjustment for your neck and back.



Power Lumbar

Individually adjustable lumbar support, which can assist in maintaining a proper sitting posture; helping to relieve tension on the spine to prevent pain or discomfort due to poor posture, arthritis or existing back problems.



Power Backrest and Power Footrest

Independently adjust the backrest and the footrest for the position that best fits the clinical needs of the client.

THE TWO STAGE LIFT MECHANISM HAS A MORE VERTICAL LIFT WITH A SLIGHT TILT COMPARED TO THE FORWARD TILT YOU FIND ON GENERIC CHAIRS.



Easy to use hand control

- Soft grip buttons
- 6 pre-defined positions
- Choose “full recline” or “full stand” with the touch of a button.

General Conditions that the Power Headrest and Lumbar adjustment can assist with:

- Kyphosis
- Lordosis
- Arthritis



PICASSO LIFT CHAIR

Customer Feedback

“I came over to your shop yesterday to try out the latest chair to hit Australia and wow what a chair! There is nothing in your shop the comes close and I’ve tried them all. Having a fused neck, arthritis and just had two lower back surgeries I am needing support while sitting. I could never get a chair to support my lower back and neck until now because of the electric adjustable lumbar support and electric adjustable head rest. Gone are the days of trying to use pillows to get the support which could never be achieved. Even the way the chair lifts you has improved as you don’t feel like you are going to land on your face.”

Therapist Feedback

“I feel this would have a good application for kyphosed clients and offer adjustability not in the typical lift chairs ranges on offer. The contoured seat base is also a good feature, to assist with lower limb and pelvic positioning, acting as a lateral thigh guide to a degree. Although often client with a higher or more complex seating needs may not be the best candidate for the standard lift chairs ranges, these will be welcome additions to provide low to moderate seating support and correction which I think will be beneficial.”

Improving your life everyday™

Sales, administration
and general enquiries
National 1300 783 783

Your Customer Service Team
customerservice@kcare.com.au

kcare.com.au