



# REHAB SHOWER COMMODOES

Improving your  
life everyday™

**KCare**  
Healthcare Solutions

# ATTENDANT PROPELLED

## REHAB SHOWER COMMUNE



\*Attendant Propelled Commode model shown with optional extras - headrest, seat, laterals, knee supports, armrest pads, neoprene footplate covers. Options sold separately.



Lateral thoracic support adjustable height mounting plates std on all models



## Features

- 400mm, 445mm, and 500mm seat widths
- Lateral thoracic support adjustable height mounting plates std on all models 240mm range (20mm increments)
- Height adjustability of the frame
- Knee and Hip guide mounting plates standard on all models 40mm horizontal mounting range (20mm increments)
- Height adjustable flip up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Standard left and right side opening on seat frame
- Seat depth adjustment
- Durable and long lasting with a stainless steel frame
- Swing away footrests
- 4 x 125mm locking swivel castors

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Junior Attendant Propelled Rehab Commode <b>KA114SJ-REHAB</b>	400mm	435mm – 465mm	360mm	500mm	905mm	575mm – 475mm	552mm – 452mm	200kg
Adult Attendant Propelled Rehab Commode <b>KA114S-REHAB</b>	445mm	435mm – 465mm	360mm	540mm	905mm	575mm – 475mm	552mm – 452mm	200kg
Adult Attendant Propelled Rehab Commode <b>KA114S50-REHAB</b>	500mm	435mm – 465mm	360mm	595mm	905mm	575mm – 475mm	552mm – 452mm	200kg

Conforms to AS/NZS 3973:2009 for safety and peace of mind

# SELF PROPELLED

## REHAB SHOWER COMMUNE



\*Self Propelled Commode model shown with optional extras - headrest, seat, laterals, knee supports, armrest pads, rotating safety arms. Options sold separately.



2 x Front 125mm  
locking swivel castors



## Features

- 400mm, 445mm, and 500mm seat widths
- Lateral thoracic support adjustable height mounting plates std on all models 240mm range (20mm increments)
- Height adjustability of the frame
- Knee and Hip guide mounting plates standard on all models 40mm horizontal mounting range (20mm increments)
- Height adjustable flip up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Standard left and right side opening on seat frame
- Seat depth adjustment
- Durable and long lasting with a stainless steel frame
- Swing away footrests
- Large self-propelling rear wheels are flush with seat top for unrestricted side transfer

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Junior Self Propelled Rehab Commode <b>KA120SRJ-REHAB</b>	400mm	435mm – 465mm	360mm	585mm	1050mm	550mm – 475mm	527mm – 452mm	200kg
Adult Self Propelled Rehab Commode <b>KA120SR-REHAB</b>	445mm	435mm – 465mm	360mm	630mm	1050mm	550mm – 475mm	527mm – 452mm	200kg
Adult Self Propelled Rehab Commode <b>KA120SR50-REHAB</b>	500mm	435mm – 465mm	360mm	685mm	1050mm	550mm – 475mm	527mm – 452mm	200kg

Conforms to AS/NZS 3973:2009 for safety and peace of mind

# PAEDIATRIC ATTENDANT PROPELLED

## REHAB SHOWER COMMODE



### Features

- 360mm and 320mm seat widths
- Lateral thoracic support adjustable height mounting plates std on all models 240mm range (20mm increments)
- Height adjustability of the frame
- Knee and Hip guide mounting plates standard on all models 40mm horizontal mounting range (20mm increments)
- Height adjustable flip up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Option of left and right side opening on seat frame (not standard on this model)
- Seat depth adjustment
- Durable and long lasting with a stainless steel frame
- Swing away footrests
- 4 x 125mm locking swivel castors

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Child Attendant Propelled Rehab Commode <b>KA114SK-REHAB</b>	360mm	320mm	360mm	485mm	750mm	575mm – 475mm	552mm – 452mm	120kg
Paediatric Attendant Propelled Rehab Commode <b>KA114SC-REHAB</b>	320mm	320mm	360mm	485mm	750mm	575mm – 475mm	552mm – 452mm	120kg

Conforms to AS/NZS 3973:2009 for safety and peace of mind

\*Paediatric Attendant Propelled Commode model shown with optional extras – lateral, seat, knee support, armrest pads, rotating safety arms, neoprene footplate covers. Options sold separately.

# PAEDIATRIC SELF PROPELLED

## REHAB SHOWER COMMODE

### Features

- 360mm and 320mm seat widths
- Lateral thoracic support adjustable height mounting plates std on all models 240mm range (20mm increments)
- Height adjustability of the frame
- Knee and Hip guide mounting plates standard on all models 40mm horizontal mounting range (20mm increments)
- Height adjustable flip up arms for unrestricted side transfer
- Adjustable tension backrest upholstery
- Option of left and right side opening on seat frame (not standard on this model)
- Seat depth adjustment
- Durable and long lasting with a stainless steel frame
- Swing away footrests
- Large self-propelling rear wheels are flush with seat top for unrestricted side transfer
- 2 x Front 125mm locking swivel castors

Description	Seat Width	Seat Depth	Backrest Height	Overall Width	Overall Depth	Floor to Seat Height Front Range	Floor to Seat Height Rear Range	Max User Weight
Child Self Propelled Rehab Commode <b>KA120SRK-REHAB</b>	360mm	320mm	360mm	630mm	900mm	550mm – 475mm	527mm – 452mm	120kg
Paediatric Self Propelled Rehab Commode <b>KA120SRC-REHAB</b>	330mm	320mm	360mm	630mm	900mm	550mm – 475mm	527mm – 452mm	120kg

Conforms to AS/NZS 3973:2009 for safety and peace of mind

# OPTIONS AVAILABLE

---

- Lateral thoracic supports
  - Headrest
  - Knee supports/ Hip Guides without adjustable link
  - Knee supports / Hip Guides with adjustable link
  - Pommel/splash guard
  - Adjustable angle footplates
  - Strapped back
  - Fixed curved padded back
  - Armrest pads
  - Rotating Safety Arms
  - Gutter arms
  - Neoprene footplate covers
  - Vinyl padded individual footplates
  - Footrest extensions
- 

**Improving your  
life everyday™**

**Sales, administration  
and general enquiries**  
National 1300 783 783

**Your Customer Service Team**  
customerservice@kcare .com.au

**kcare.com.au**