



# Etac Cloo toilet seat raiser

Designed to assist users and carers with every day routines.

 **etac**<sup>®</sup>

Creating Possibilities

# Etac Cloo toilet seat raiser

Cloo is a height adjustable toilet seat raiser which combines flexibility and simplicity. It has a discreet and soft design to blend in with your bathroom interior.

## Adjustable angle and height

Cloo can be adjusted in three different heights 6 cm, 10 cm or 14 cm. It can also be set with an anterior tilt which facilitates sitting down and standing up for users with reduced mobility in hip or knee joints.

## Functional arm supports

Cloo is available with or without detachable arm supports. The arm supports are foldable but can be fixed in a folded down position if required. Cloo can be used with only one arm support if space around the toilet is limited.

## Easy to clean and mount

To aid cleaning Cloo has soft plain surfaces and the seat stays folded in the upright position. The seat can also be easily detached should this be preferred. Mounting is quick and tool-less.



### With or without arm supports

Etac Cloo is available with or without arm supports.



### Interior tilt

Etac Cloo can be set with an anterior tilt which facilitates sitting down and standing up for users with reduced mobility in hip or knee joints.



### Easy to clean

The seat can be easily detached for easy cleaning.



## Accessory



**Description** Arm supports, pair  
**Item no.** 80301222

**Material** Polypropylene, powder coated steel tubing



**Description** Heightening spacers, 14 cm 5½"  
**Item no.** 80301221

Standard height adjustability: 6 and 10 cm 2¼" and 4" as standard. These are needed if 14 cm 5½" is required.

**Material** Polyamide

Description	Dimensions	Weight	Item no.
Cloo Toilet Seat Raiser	Seat width 40 cm 15¾" Total width 60 cm 23½"	3.3 kg 7.28 lbs	80301204 ✓
Cloo Toilet Seat Raiser with Armrests	Seat width 40 cm 15¾" Total width 60 cm 23½" Height of Armrests from seat 21 cm 8¼" Width between Armrests 49 cm 19¼"	4.8 kg 10.58 lbs	80301214 ✓

### Material

Seat: Polypropylene  
Attachment: Aluminium

Clean the product with a solvent-free cleaning agent (pH level 6–9), or with a 70% disinfectant solution. Can be decontaminated at max. 85°C / 185°F if so required.



max 85°C max 185°F  
for maximum of 3 minutes


max 150 kg max 330 lbs





 **etac**<sup>®</sup>  
Creating Possibilities

“Everyone should be able to live a free and independent life and pursue their dreams regardless of any physical circumstances.”

 = Standard Stock Item